

WEGO



WEIGAO MEDICAL INTERNATIONAL CO., LTD.

ADD: No.1 Weigao Road,Hi-Tech Zone,Weihai,Shandong,China
TEL: +86 631 5660699 FAX: +86 631 5622563 P.C.: 264209
E-mail: info@weigaogroup.com exp@weigaogroup.com
Http: //www.wegomedintl.com/

This copyright belongs to Wego Medical International Co., Ltd. Our company reserves the right to investigate legal liability in case of infringement.

2022.9 EDITION



Enteral Nutrition Feeding Set

>>> ENTERAL - SYRINGES

Syringes from 1 ml to 50/60 ml for the enteral administration of medication, flushing of ducts, gastric residuals and breast milk or for enteral administration by gravity, bolus, or syringe pump.



>>> Purple plunger, double markings for 1 ml, 3ml, 5ml, 10ml, 20ml, 30ml, and 50/60ml.

Enteral - Jeringas				
CODE NUMBER	VOLUMEN (ML)	TIP POSITION	DISPENSOR UNITS	UNIT PER BOX
WG01-010-001-001	1	Centered	200	3000
WG01-010-001-002	3	Centered	200	1800
WG01-010-001-003	5	Centered	100	1600
WG01-010-001-004	10	Centered	100	1200
WG01-010-001-005	20	Centered	100	600
WG01-010-001-006	30	Centered	50	400
WG01-010-001-007	50/60	Centered	20	200

>>> Enteral meets the latest international safety requirements defined by the ISO Committee.



>>> ENFIT COMPONENTS



WG01-010-002-001
Female connector, cap on strap



WG01-010-002-002
Y-connector, Male sideport with cap on strap



WG01-010-002-003
Non-vented Enplus cross spike with vented cap



WG01-010-002-004
Transition adapter for tube



WG01-010-002-005
Two-way tube+magnetic connector



WG01-010-002-006
Fluid exit connector

ENTERAL FEEDING BAG Twins Bag



Enteral Feeding Sets

Product name

Enteral feeding bag

For pump use

Code: WG01-010-003

Speification parameter

NO	VOLUME(ML)
WG01-010-003-001	500
WG01-010-003-002	600
WG01-010-003-003	1000
WG01-010-003-004	1200
WG01-010-003-005	1500

Product introduction

- PVC material, used with infusion pumps
- Tube length 130 cm, tube diameter 4.3 mm. Customizable bag body and tube size.
- DEHP-free available
- EO sterilized.

ENTERAL FEEDING BAG

Gravity using



Enteral Feeding Sets

Product name

Enteral feeding bag

For gravity feeding

Code: WG01-010-004

Speification parameter

NO	VOLUME(ML)
WG01-010-004-001	500
WG01-010-004-002	600
WG01-010-004-003	1000
WG01-010-004-004	1200
WG01-010-004-005	1500

Product introduction

- Special PVC material for gravity feeding device.
- Tube length 130 cm, tube diameter 4.3 mm. Customizable bag body and tube size.
- DEHP-free available
- EO sterilized.

ENTERAL FEEDING SET PUNCTURE SPIKE/CAP CONNECTOR



Enteral Feeding Sets

Product name

Enteral feeding set


Puncture Spike/Cap Connector

Code: WG01-010-005

Product introduction

- Product: Disposable enteral feeding set(Gravity feeding)
- Material: Nontoxic medical grade PVC.
- Package: Individual paper-poly bag
- DEHP-free available

ENTERAL FEEDING SET PUNCTURE SPIKE/CAP CONNECTOR



Enteral Feeding Sets

Product name

Enteral feeding set
Puncture Spike/Cap Connector
Code: WG01-010-006

Product introduction

- Product: Disposable enteral feeding set(Pump feeding)
- Applicable: Peristaltic pump/Rotatory pump
- No-toxic medical grade PVC and silicon.
- Package: Individual paper-poly bag
- DEHP-free available