

WEGO



WEIGAO MEDICAL INTERNATIONAL CO., LTD.

ADD: No.1 Weigao Road,Hi-Tech Zone, Weihai, Shandong, China

TEL: +86 631 5660699 FAX: +86 631 5622563 P.C.: 264209

E-mail: info@weigaogroup.com exp@weigaogroup.com

Http: [//www.wegomedintl.com/](http://www.wegomedintl.com/)

This copyright belongs to Weigao Medical International Co., Ltd. Our company reserves the right to investigate legal liability in case of infringement.

2023.2 EDITION



Infusion Sets with Precision Filters for Single Use

About WEGO Group

Established in 1988, with more than 30 years of dedication to medical devices and pharmaceuticals production, Wego continuously strives for excellence in developing and delivering quality products and services in healthcare. We focus on Clinical Nursing, Orthopedics, Blood Purification, IVD Supply, Wound Care, Surgery Consumables, Medical Imaging, Intervention, Medical Equipment, and Service. Based on the value of "Your health We care," we tirelessly improve our products, which are recognized and exported to over 150 countries and regions. WEGO's products are FDA-approval and are certified to CE, ISO, GMP, and IVDD.

Today, Wego's 30,000 employees across the world set our mission, making the world a better place by bringing life-changing health technologies and innovative medical system solutions to the people.





WG01-001-001	ILLUSTRATION
SM-001	
SM-002	
SM-003	

Chamber Infusion Sets (WG01-001-001)

Air-vented spike

- For Infusion bags With solution filter
- 20 Drops= 1 ± 0.1 ml
- Flexible Drip Chamber(6ml)
- Roller Clamp
- Slip or Lock Adapter
- Y Injection Site (Latex or Latex free)
- For Gravity Use
- without needle
- Sterile




WG01-001-005	ILLUSTRATION
SM-013	
SM-014	
SM-015	

Chamber Infusion Sets (WG01-001-005)

Air-vented spike

- For Infusion bags With solution filter
- 20 Drops= 1 ± 0.1 ml
- Flexible Drip Chamber(6ml)
- Roller Clamp
- Slip or Lock Adapter
- Latex Injection Site With/ without needle
- For Gravity Use
- Sterile



WG01-001-002	ILLUSTRATION
SM-004	
SM-005	
SM-006	

Chamber Infusion Sets (WG01-001-002)

Air-vented spike

- For Infusion bags and bottle With solution filter
- 20 Drops= 1 ± 0.1 ml
- Flexible Drip Chamber(6ml)
- Roller Clamp
- Slip or Lock Adapter
- Y Injection Site (Latex or Latex free)
- For Gravity Use
- With needle
- Sterile



WG01-001-006	ILLUSTRATION
SM-016	
SM-017	
SM-018	

Chamber Infusion Sets (WG01-001-006)

Air-vented spike

- For Infusion bags With solution filter
- 20 Drops= 1 ± 0.1 ml
- Flexible Drip Chamber(6ml)
- Roller Clamp
- Slip or Lock Adapter
- Latex/Latex-free YInjection Site With/ without needle
- For Gravity Use
- Sterile



WG01-001-003	ILLUSTRATION
SM-007	
SM-008	
SM-009	

Chamber Infusion Sets (WG01-001-003)

Air-vented spike

- For Infusion bags With solution filter
- 20 Drops= 1 ± 0.1 ml
- Flexible Drip Chamber(6ml)
- Roller Clamp
- Slip or Lock Adapter
- Latex/Latex-free Y Injection Site With/ without needle
- For Gravity Use
- Sterile



WG01-001-007	ILLUSTRATION
BI-001	
BI-002	
BI-003	

Chamber Infusion Sets (WG01-001-007)

With air-vent spike

- For Infusion bags/bottle
- 20 Drops= 1 ± 0.1 ml With solution filter
- Flexible Drip Chamber(8-10ml)
- Precision regulator
- Slip or Lock Adapter
- Latex/Latex-free YInjection Site With/ without needle
- For Gravity Use
- Sterile



WG01-001-004	ILLUSTRATION
SM-010	
SM-011	
SM-012	

Chamber Infusion Sets light resistant (WG01-001-004)

Air-vented spike

- For Infusion bags With solution filter
- 20 Drops= 1 ± 0.1 ml
- Flexible Drip Chamber(6ml)
- Roller Clamp
- Slip or Lock Adapter
- Latex/Latex-free Y Injection Site With/ without needle
- For Gravity Use
- Sterile



WG01-001-008	ILLUSTRATION
BI-004	
BI-005	
BI-006	

Chamber Infusion Sets Light Resistant (WG01-001-008)

Air-vented spike Light-proof type

- For Infusion bags With solution filter
- 20 Drops= 1 ± 0.1 ml
- Flexible Drip Chamber(8-10ml)
- ABSRoller Clamp
- Slip or Lock Adapter
- Latex/Latex-free YInjection Site With/ without needle
- For Gravity Use
- Sterile



WG01-001-009	ILLUSTRATION
BI-007	
BI-008	
BI-009	

Chamber Infusion Sets (WG01-001-009)

Air-vented spike

- For Infusion bags With solution filter
- 60 Drops= 1 ± 0.1 ml
- Flexible Drip Chamber(8-10ml)
- With Clip
- Slip or Lock Adapter
- Latex/Latex-free YInjection Site With/ without needle
- For Gravity Use
- Sterile




WG01-001-010	ILLUSTRATION
BI-010	
BI-011	
BI-012	

Chamber Infusion Sets (WG01-001-010)

Air-vented spike

- For Infusion bags With solution filter
- 20 Drops= 1 ± 0.1 ml
- Flexible Drip Chamber(8-10ml)
- Roller Clamp
- Slip or Lock Adapter
- Protecting type, Latex/Latex-free YInjection Site With/ without needle
- For Gravity Use
- Sterile



WG01-001-011	ILLUSTRATION
IN-001	
IN-002	
IN-003	

Infusion Sets (WG01-001-011)

Without air-vent spike

- For Infusion bags With solution filter
- 20 Drops= 1 ± 0.1 ml
- Flexible Drip Chamber
- Roller Clamp
- Slip or Lock Adapter
- For Gravity Use With/ without needle
- For Gravity Use
- Sterile



WG01-001-012	ILLUSTRATION
IN-004	
IN-005	
IN-006	

Infusion Sets Light Resistant (WG01-001-012)

- For Infusion bags
- 20 Drops= 1 ± 0.1 ml
- Flexible Drip Chamber
- Roller Clamp
- Slip or Lock Adapter
- Latex/Latex-free Injection Site
- For Gravity Use
- Sterile



WG01-001-013	ILLUSTRATION
DO-001	
DO-002	
DO-003	

Double Chamber Infusion Sets (WG01-001-013)

- For Infusion bags With air-vent With solution filter
- 20 Drops= 1 ± 0.1 ml
- Flexible Drip Chamber
- SafetyRoller Clamp
- Rotated Lock Adapter
- Positive pressure needle-free connector
- For Gravity Use
- with/without needle
- Sterile



WG01-001-014	ILLUSTRATION
DO-004	
DO-005	
DO-006	

Double Chamber Infusion Sets (WG01-001-014)

- For Infusion bags/bottle With air-vent With solution filter
- 20 Drops= 1 ± 0.1 ml
- Flexible Drip Chamber
- SafetyRoller Clamp
- Rotated Lock Adapter
- Latex/Latex-free YInjection Site With/ without needle
- For Gravity Use
- Sterile



WG01-001-015	ILLUSTRATION
SO-001	
SO-002	
SO-003	

Soft Burette(WG01-001-015)




- 100 and 150ml Graduated Burette With Shut-off Valve
- 60 Drops= 1 ± 0.1 ml
- Flexible Drip Chamber with Filter
- Roller Clamp
- Slip or Lock Adapter with Protective Cap
- Latex/Late-free Y Injection Site With/without needle
- For Gravity Use
- Sterile



WG01-001-016	ILLUSTRATION
IN-007	

Infusion Sets (WG01-001-016)

- For Infusion bags
- 20 Drops= 1 ± 0.1 ml
- Flexible Drip Chamber
- Roller Clamp
- Slip or Lock Adapter
- Latex/Latex-free Injection Site
- For Gravity Use
- Sterile



Extension Tube

- 36" (91 cm), Appx 3.8mL
- Microbore ext set with clamp, fixed male luer, non-DEHP
- 0.8mm x 2.4mm
- Sterile



Two Way Extension Tube


- 7" (18 cm), Appx 0.3mL
- Smallbore, fixed male luer lock, clamp, non-DEHP
- 0.8mm x 2.4mm
- Sterile



- 5" (12 cm), Appx 0.3 mL
- Microbore, t-connector, clamp, non-DEHP
- 0.8mm x 2.4mm
- Sterile



FR62-3 FR62-3
Flow Regulator



Water Blocking Protective cap

Y-injection Port



Flow Regulator



Three-way Valve

