

WEGO



WEIGAO MEDICAL INTERNATIONAL CO., LTD.

ADD: No.1 Weigao Road,Hi-Tech Zone, Weihai,Shandong,China
TEL: +86 631 5660699 FAX: +86 631 5622563 P.C.: 264209
E-mail: info@weigaogroup.com exp@weigaogroup.com
Http: //www.wegomedintl.com/

This copyright belongs to Wego Medical International Co., Ltd. Our company reserves the right to investigate legal liability in case of infringement.

2022.11 EDITION



Disposable Hypodermic Syringe and Needle

About WEGO Group

Established in 1988, with more than 30 years of dedication to medical devices and pharmaceuticals production, Wego continuously strives for excellence in developing and delivering quality products and services in healthcare. We focus on Clinical Nursing, Orthopedics, Blood Purification, IVD Supply, Wound Care, Surgery Consumables, Medical Imaging, Intervention, Medical Equipment, and Service. Based on the value of "Your health We care," we tirelessly improve our products, which are recognized and exported to over 150 countries and regions. WEGO's products are FDA-approval and are certified to CE, ISO, GMP, and IVDD.


Today, Wego's 30,000 employees across the world set our mission, making the world a better place by bringing life-changing health technologies and innovative medical system solutions to the people.



Disposable Syringe (WG01-003-001)

- Material: Polypropylene
- Parts: Barrel+Piston+Plunger+Needle(or w/o needle) Luer Lock or Luer Slip
- Packing: Blister/PE
- Piston materials: Latex, Polyisoprene, Butyl
- Sterilization: EO

Code Number	Volume (ml/cc)	Needle	Luer Lock	Luer Slip	Packing
WGS1WL	1	W	L	S	Blister/ PE
WGS1WOS		W/O	L	S	Blister/ PE
WGS1WS		W	L	S	Blister/ PE
WGS1WOL	2	W	L	S	Blister/ PE
WGS2WL		W/O	L	S	Blister/ PE
WGS2WOS		W	L	S	Blister/ PE
WGS2WS	3	W	L	S	Blister/ PE
WGS2WOL		W/O	L	S	Blister/ PE
WGS3WL		W	L	S	Blister/ PE
WGS3WOS	5	W	L	S	Blister/ PE
WGS3WS		W/O	L	S	Blister/ PE
WGS3WOL		W	L	S	Blister/ PE
WGS5WL	10	W	L	S	Blister/ PE
WGS5WOS		W/O	L	S	Blister/ PE
WGS5WS		W	L	S	Blister/ PE
WGS5WOL	20	W	L	S	Blister/ PE
WGS10WL		W/O	L	S	Blister/ PE
WGS10WOS		W	L	S	Blister/ PE
WGS10WS	30	W	L	S	Blister/ PE
WGS10WOL		W/O	L	S	Blister/ PE
WGS20WL		W	L	S	Blister/ PE
WGS20WOS	50	W	L	S	Blister/ PE
WGS20WS		W/O	L	S	Blister/ PE
WGS20WOL		W	L	S	Blister/ PE
WGS30WL	60	W	L	S	Blister/ PE
WGS30WOS		W/O	L	S	Blister/ PE
WGS30WS		W	L	S	Blister/ PE
WGS30WOL	100	W	L	S	Blister/ PE
WGS50WL		W/O	L	S	Blister/ PE
WGS50WOS		W	L	S	Blister/ PE
WGS50WS	60	W	L	S	Blister/ PE
WGS50WOL		W/O	L	S	Blister/ PE
WGS60WL		W	L	S	Blister/ PE
WGS60WOS	100	W	L	S	Blister/ PE
WGS60WS		W/O	L	S	Blister/ PE
WGS60WOL		W	L	S	Blister/ PE
WGS100WL	100	W	L	S	Blister/ PE
WGS100WOS		W/O	L	S	Blister/ PE
WGS100WS		W	L	S	Blister/ PE
WGS100WOL		W/O	L	S	Blister/ PE




Disposable Insulin Syringe (WG01-003-002)

- Material: Polypropylene
- Parts: Barrel+Piston+Plunger+Fixed Needle+Cap
- Packing: Individual Blister/Self-Sealed Package
- Piston materials: Latex, Polyisoprene, Butyl
- Sterilization: EO

Code Number	Volume (ml/cc)	Needle Size		Packing	Colour Identification
		OD	L		
WGIN5J1	1	29G	1/2	Blister/ PE	Orange
WGIN5J5	0.5	30G	3/8		
WGIN5J3	0.3	31G	5/16		
WGIN5J1	1	29G	1/2	Blister/ PE	Red
WGIN5J5	0.5	30G	3/8		
WGIN5J3	0.3	31G	5/16		

* Customized needle size is available



Disposable Perfusor Syringe (WG01-003-003)

- Material: Polypropylene
- Parts: Barrel+Piston+Plunger
- Packing: Blister/PE
- Piston materials: Latex, Polyisoprene, Butyl
- Sterilization: EO
- Type: Conventional/Light-Resistant
- Pressure: 200pa
- Size: 50/60 ML(cc)

Pre-filled Flush Syringe (WG01-003-004)

- Sterilization method: high temperature moist heat sterilization
- Main materials: polypropylene, butyl rubber

Specification	Single PE bag	Box size	Capacity of a box	Carton size	G.W
3ml	30pcs/box	26.6*19*9.1cm	480pcs/carton	530*385*390mm	8kg
5ml	30pcs/box	22.5*16*9.5cm	480pcs/carton	475*345*415mm	10kg
10ml	30pcs/box	22*22*9.5cm	480pcs/carton	475*465*415mm	14kg



Disposable Sterile Rinsers with Catheter Tip (WG01-003-005)

- Material: Polypropylene
- Parts: Barrel+Piston+Plunger+Cap
- Packing: Blister/PE
- Piston materials: Latex, Polyisoprene, Butyl
- Sterilization: EO
- Size: 50/60 ML(cc)



Feeding Syringe (WG01-003-006)

Feeding syringe are designed specifically for enteral processes. These processes include initial tube placement, flushing, irrigation and more. The ENFit connector reduces the risk of misconnections to the tubing. Also, the body is clear for easy measuring against the clearly marked graduated length markings.


- Material: Polypropylene
- Parts: Barrel+Piston+Plunger+Cap
- Packing: Blister/PE
- Piston materials: Latex, Polyisoprene, Butyl
- Sterilization: EO
- Size: 3/5/ 10/20/30/50 ML(cc)



Disposable Pump Syringe (WG01-003-007)

- Material: Polypropylene
- Parts: Barrel+Piston+Plunger
- Packing: Blister/PE
- Piston materials: Latex, Polyisoprene, Butyl
- Sterilization: EO

Code Number	Volume (ml/cc)	Needle		Luer Lock	Packing
WGPS3	3	W	W/O	L	Blister/ PE
WGPS5	5	W	W/O	L	Blister/ PE
WGPS10	10	W	W/O	L	Blister/ PE
WGPS20	20	W	W/O	L	Blister/ PE
WGPS30	30	W	W/O	L	Blister/ PE
WGPS50	50	W	W/O	L	Blister/ PE
WGPS60	60	W	W/O	L	Blister/ PE



Disposable Hypodermic Needle



Code Number	Specifications		Colour		Regular Needle Length							
					13mm	16mm	19mm	25mm	32mm	38mm	50mm	
	OD(*)	GAUGE	1/2"(A)	5/8"(B)	3/4"(C)	1"(D)	1 1/4"(E)	1 1/2"(F)	2"(G)			
WGZ30	0.30	30G	Light yellow									
WGZ29	0.33	29G	Red									
WGZ28	0.36	28G	Turquoise									
WGZ27	0.40	27G	Grey									
WGZ26	0.45	26G	Brown									
WGZ25	0.50	25G	Orange									
WGZ24	0.55	24G	Medium purple									
WGZ23	0.60	23G	Sky blue									
WGZ22	0.70	22G	Black									
WGZ21	0.80	21G	Dark green									
WGZ20	0.90	20G	Yellow									
WGZ19	1.10	19G	Cream yellow									
WGZ18	1.20	18G	Pink									
WGZ17	1.5	17G	Purple									
WGZ16	1.6	16G	White					White		White		

* Customized needle size is available

Disposable Light-resistant Syringe (WG01-003-008)

- Material: Polypropylene
- Parts: Barrel+Piston+Plunger
- Packing: Blister/PE
- Piston materials: Latex, Polyisoprene, Butyl
- Sterilization: EO


Code Number	Volume (ml/cc)	Needle		Luer Lock	Packing
WGLRS3	3	W	W/O	L	Blister/ PE
WGLRS5	5	W	W/O	L	Blister/ PE
WGLRS10	10	W	W/O	L	Blister/ PE
WGLRS20	20	W	W/O	L	Blister/ PE
WGLRS30	30	W	W/O	L	Blister/ PE
WGLRS50	50	W	W/O	L	Blister/ PE



Disposable Fixed Dose Vaccine Syringe (WG01-003-009)

- Material: Polypropylene
- Parts: Barrel+Piston+Plunger+Fastener+Fixed Needle+Cap
- Packing: Blister/PE
- Piston materials: Latex, Polyisoprene, Butyl
- Sterilization: EO

Code Number	Volume (ml/cc)	Needle Size		Packing
		OD	L	
WGV501	1	26G	1/2"	Blister/ PE
WGV502	0.9			Blister/ PE
WGV503	0.8			Blister/ PE
WGV504	0.7	25G	5/8"	Blister/ PE
WGV505	0.6			Blister/ PE
WGV506	0.5	24G	5/8"	Blister/ PE
WGV507	0.4			Blister/ PE
WGV508	0.3	23G	1"	Blister/ PE
WGV509	0.2			Blister/ PE
WGV51	0.1			Blister/ PE



Disposable Safety Hypodermic Needle

- Material: Polypropylene+Stainless Steel
- Parts: Needle+Cap
- Packing: Blister
- Sterilization: EO
- Needle Size: 18G-27G

