

# Disposable Oral Enteral Syringe

ENFit compliant

Sterilization

DEHP/Latex Free



# Product Specification

From 1ml to 60ml, meet the requirement of different sizes



## Applicable scene

Safe for neonates, pets and adult patients



# Product Composition

**Cap**

Prevent syringes from leaking

**Barrel**

Transparent barrel easy to observe the liquid level and bubbles

**PISTON**

Leakage-proof plunger guarantee, special design to prevent the plunger from slipping out

**PLUNGER**

The Purple-colored plunger rod is designed to help easily differentiate between a traditional parenteral syringe and one intended for use with enteral feeding



# Low dose tip design

Improve dose accuracy (1ml and 2.5ml)



# Clear Scale

A clear barrel lets you easily view the medication



# Syringe Accessories

WEGO has put together a large range of accessories to help with the delivery of feeds and medicines using the ENFIT syringe feeding system



**WG01-010-002-001**  
Female connector, cap on strap



**WG01-010-002-002**  
Y-connector, Male sideport with cap on strap



**WG01-010-002-003**  
Non-vented Enplus cross spike with vented cap



**WG01-010-002-004**  
Transition adapter for tube



**WG01-010-002-005**  
Two-way tube+magnetic connector



**WG01-010-002-006**  
Fluid exit connector