

HYPODERMIC NEEDLE



**INDIVIDUAL
PACKAGE**



**EO
STERILIZATION**



**GENERAL
INTERFACE**

GENERAL INTERFACE

Adaptable to different interfaces and different types of syringes, easy to operate



LUER LOCK



LUER SLIP

PRODUCT DISPLAY



NEEDLE SHARP

Fast needle puncture,
less pain and less tissue damage



Three blade structure

High quality stainless steel

INDIVIDUAL PACKAGE

Soft and hard blister packing



Hypodermic Needle



SIZE		COLOR		LENGTH OF REGULAR NEEDLES					
O.D (mm)	GUAGE	COLOR CODE		13mm 1/2"	16mm 5/8"	19mm 3/4"	25mm 1"	32mm 1 1/4"	38mm 1 1/2"
1.60	16G		White						
1.20	18G		Pink						
1.00	19G		Beige						
0.90	20G		Yellow						
0.80	21G		Green						
0.70	22G		Black						
0.60	23G		Blue						
0.55	24G		Purple						
0.50	25G		Orange						
0.45	26G		Brown						
0.40	27G		Grey						
0.36	28G		Blue-Green						
0.33	29G		Red						
0.30	30G		Light Yellow						

*Other specifications are available according to customer requests.