

# INFUSION SETS

Mainly used for the transportation of liquid medicine for injured patients



EU CE CERTIFICATION



ISO CERTIFICATION



FDA 510K



## SHARP NEEDLE

Stainless steel needle (With or Without), hypodermic needle and butterfly type needle are available.

---



## TRANSPARENT TUBE

Soft, kink resistant medical grade PVC tubing with length 150 cm or customized.

---



## FLOW REGULATOR

Accurately control the flow rate, normal shape and horned shape are available.

---





## CONNECTOR CONE

Made of PE or PVC, luer lock and luer slip are available.

## DRIP CHAMBER

Made of medical grade PVC, with or without filter. 20 drops/ml and 60 drops/ml are available.



LIQUID FILTER



## RELATED ACCESSORIES



# PRODUCT SPECIFICATIONS

## INFUSION SETS FOR SINGLE USE (WITH NEEDLE)

---

[Learn more >](#)



## INFUSION SETS FOR SINGLE USE (WITHOUT NEEDLE)

---

[Learn more >](#)



## INFUSION SETS WITH PRECISION FILTERS FOR SINGLE USE (WITH NEEDLE)

---

[Learn more >](#)



## INFUSION SETS WITH FILTERS FOR SINGLE USE (WITHOUT NEEDLE)

---

[Learn more >](#)



## INFUSION SETS WITH BURETTE FOR SINGLE USE (WITHOUT NEEDLE)

---

[Learn more >](#)



## LIGHT-RESISTANT INFUSION SETS FOR SINGLE USE (WITHOUT NEEDLE)

---

[Learn more >](#)

