

SYRINGE SERIES


MULTIPLE TYPES

OEM&ODM


COMPLETE CERTIFICATE: CE、ISO9001、ISO13485



/ INSULIN SYRINGE



/ IRRIGATION SYRINGE



/ ENFIT SYRINGE



/ STERILE SYRINGE



/ LIGHT RESISTANT SYRINGE



/ PUMP SYRINGE



/ VACCINE SYRINGE

PRODUCT COMPOSITION

SCIENTIFIC DESIGN

HYPODERMIC NEEDLE

MANY SIZES OF NEEDLES ARE AVAILABLE

PISTON

DOUBLE-LAYER RUBBER STOPPER HAS GOOD AIR TIGHTNESS

TRANSPARENT COAT

TRANSPARENT COAT EASY TO OBSERVE THE LIQUID LEVEL AND BUBBLES

PLUNGER

ANTI-SLIP STRUCTURE TO PREVENT THE CORE ROD FROM ACCIDENTALLY SLIPPING OFF





TWO SPECIFICATIONS

LUER LOCK/LUER SLIP

LUER LOCK




LUER SLIP



CLEAR SCALE

STRONG INK ADHESION







GOOD SEALING

THE RUBBER STOPPER HAS A HIGH DEGREE OF ADHESION TO THE TUBE



COMPLETE SPECIFICATIONS


MULTIPLE SPECIFICATIONS, MORE CHOICES





SHOW DETAILS

DESIGN EVERY DETAIL WITH CARE




STERILE INDIVIDUAL PACKAGING

EO STERILIZATION

BLISTER



PE BAG



CERTIFICATE DISPLAY

QUALITY IS TRUSTWORTHY

