

WEGO

Sterilizer Series

To be the strongest in China,best-in-class globally,the most respected innovative medical devices and pharmaceutical company



WEIGAO MEDICAL INTERNATIONAL CO., LTD.

Add: NO.1 Weigao Road, Weigao Industry Zone, Weihai, Shandong Province, China.
Tel: +866315717597 FAX: +86 631 5622563
Email: shaxiaoyu@weigaogroup.com
Web: www.weigaogroup.com

This copyright belongs to Wego Medical International Co., Ltd. Our company reserves the right to investigate legal liability in case of infringement.


2023.4 EDITION

Introduction of WEGO

The WEGO was founded in 1988, a leading medical device manufacturer and has been ranked Top 500 enterprises in China, covering a manufacturing area of 6,000,000 square meters including sterilized workshops with an area of 2,600,000 square meters that comply with GMP standards.

It has 60 subsidiaries with more than 30,000 employees and currently produces more than 500 types of medical products with more than 80,000 specifications, obtaining the certifications of CE, ISO13485, China Compulsory Certification, and products registration of FDA 510K from Food and Drug Administration (FDA).

As a famous infection control solution provider, the WEGO is a leader and pioneer company in the design and manufacture of sterilizers, up-to-date design in compliance with high-quality standards and providing a safe reliable sterile processing and re-use for the items, widely used in dental, laboratory, clinic, hospital and pharmaceutical factory.





HORIZONTAL CLASS B PRESSURE STEAM STERILIZER



APPLICATION

MSG.BSD series horizontal Class B steam sterilizer is suitable for unwrapped, wrapped, solid, hollow A, hollow B, lumen instruments and implanted instrument to be sterilized, widely used in life science, OR, clinic, lab, hospital and pharmacy.

MAIN FEATURE

- Microcomputer & PLC control system, touch screen displaying
- Internal printer
- Automatic with pre-set programs
- Built-in independent rapid steam generator
- Built-in vacuum pump and vacuum drying system
- Single & Door for choice
- Sterilization temperature range: 105 C ~ 138 C
- Limit vacuum: -0.096Mpa
- Temperature uniformity: $\leq \pm 0.5$ C
- Chamber material: Stainless steel
- Accessory: Stainless steel sterilizing grille

MAIN PARAMETER

Model	Chamber volume	Chamber dimensions (Φ×L)mm	Overall dimensions (L×W×H) mm	Net weight	Power supply
MSG.BSD60L	60L	386×500	985×600×1450	105kg	7.0KW,AC380V/220V, three phase
MSG.BSD80L	80L	386×700	985×600×1450	120kg	7.0KW,AC380V/220V, three phase
MSG.BSD100L	100L	386×860	1150×600×1450	140kg	7.0KW,AC380V/220V, three phase
MSG.BSD135L	135L	500×700	1000×770×1750	225kg	14KW,AC380V/220V, three phase
MSG.BSD150L	150L	500×780	1080×770×1750	240kg	14KW,AC380V/220V, three phase
MSG.BSD185L	185L	500×950	1250×770×1750	255kg	14KW,AC380V/220V, three phase
MSG.BSD200L	200L	500×1000	1300×770×1750	265kg	14KW,AC380V/220V, three phase
MSG.BSD300L	300L	632×1000	1350×900×1850	350kg	22KW,AC380V/220V, three phase
MSG.BSD450L	450L	632×1450	1800×900×1850	450kg	32KW,AC380V/220V, three phase





HORIZONTAL CLASS N PRESSURE STEAM STERILIZER



APPLICATION

MSG.NSD series horizontal class N steam sterilizer is suitable for unwrapped, wrapped, textile, rubber, liquid, waste to be sterilized, widely used in clinic, lab and medical health center.

MAIN FEATURE

- Microcomputer & PLC control system, touch screen displaying
- Automatic with pre-set programs
- Internal printer
- Built-in independent rapid steam generator
- Single door & double door for choice
- Sterilization temperature range: 105℃ ~ 138℃
- Temperature uniformity: $\pm 0.5^{\circ}\text{C}$
- Chamber material: Stainless steel
- Accessory: Stainless steel sterilizing grille

MAIN PARAMETER

Model	Chamber volume	Chamber dimensions (Φ×L)mm	Overall dimensions (L×W×H) mm	Net weight	Power supply
MSG.NSD60L	60L	386×500	985×600×1450	90kg	6.6KW, AC380V/220V, three phase
MSG.NSD80L	80L	386×700	985×600×1450	115kg	6.6KW, AC380V/220V, three phase
MSG.NSD100L	100L	386×860	1150×600×1450	125kg	6.6KW, AC380V/220V, three phase
MSG.NSD135L	135L	500×700	1000×770×1750	210kg	12KW, AC380V/220V, three phase
MSG.NSD150L	150L	500×780	1080×770×1750	225kg	12KW, AC380V/220V, three phase
MSG.NSD185L	185L	500×950	1250×770×1750	240kg	12KW, AC380V/220V, three phase
MSG.NSD200L	200L	500×1000	1300×770×1750	250kg	12KW, AC380V/220V, three phase
MSG.NSD300L	300L	632×1000	1350×900×1850	325kg	20KW, AC380V/220V, three phase
MSG.NSD450L	450L	632×1450	1800×900×1850	420kg	30KW, AC380V/220V, three phase

