

# WEGO

LEADING BRAND OF MEDICAL  
DEVICES IN CHINA



**WEIGAO MEDICAL INTERNATIONAL CO., LTD.**

ADD: No.1 Weigao Road, Hi-Tech Zone, Weihai, Shandong, China

TEL: +86 631 5660699 FAX: +86 631 5622563 P.C.: 264210

E-mail: [info@weigaogroup.com](mailto:info@weigaogroup.com) [exp@weigaogroup.com](mailto:exp@weigaogroup.com)

Http: [//www.wegomedical.com/](http://www.wegomedical.com/)

This copyright belongs to Wego Medical International Co., Ltd. Our company reserves the right to investigate legal liability in case of infringement.

2024.2 EDITION

CE 0123

FDA





ICU Life Support







### WG01-005-001 Infusion Pump

Power supply: AC 100-240V, 50/60Hz 3.2"  
color LCD screen  
Dual CPU. Rate mode, time mode, body weight mode  
Compatible with all standard infusion sets.  
KVO and bolus function  
Various audible and visual alarms  
Optional: Wireless



### WG01-005-002 Single Channel Syringe Pump

Power supply: AC 100-240V, 50/60Hz  
compatible with 5,10,20,30,50/60 syringes of all brands  
Advanced design, anti-shock, Ingress protection IP44  
Automatically identify syringe type  
Dual CPU. Rate mode, time mode, body weight mode  
Various audible and visual alarms  
Optional: Wireless



### WG01-006 High Flow Heated Respiratory Humidifier

High performance turbine  
4.3" touch screen  
MR290 humidification tank, could reach 31°C in 10 minutes  
1000 hrs filter change reminder  
Water tank can add water automatically  
168h graphical and tabular trend  
Various audible and visual alarms



### WG01-007 Defibrillator

Manual defibrillation, AED, Pacer, 5-lead ECG, RESP, Thermal Recorder  
8.4" TFT screen  
Power supply: 100 - 240V~, 50/60Hz±1Hz  
Battery charging: Less than 2 hours to 80% and less than 3 hours to 100%  
Various audible and visual alarms  
Data storage: 2000 groups NIBP, 48 hours Wave Review  
Manual Mode, AED Mode, Synchronous defibrillation



### WG01-005-003 Double Channel Syringe Pump

Power supply: AC 100-240V, 50/60Hz  
Compatible with 5,10,20,30,50/60 syringes of all brands.  
Advanced design, anti-shock, Ingress protection IP44  
Automatically identify syringe type  
Dual CPU. Rate mode, time mode, body weight mode, online mode.  
Various audible and visual alarms  
Optional: Wireless



### WG01-008 Patient Monitor

3-lead ECG, RESP, SpO2, NIBP, HR, temp  
12.1" color TFT screen  
Power supply: 100-240V~, 50/60Hz±1Hz  
Various audible and visual alarms  
Data storage: 2000 groups NIBP, 6 hours wave review  
Built in recorder, thermal array  
Optional: Dual-IBP, EtCO2, 5/12-lead ECG