

WEGO

LEADING BRAND OF MEDICAL
DEVICES IN CHINA



WEIGAO MEDICAL INTERNATIONAL CO., LTD.

ADD: No.1 Weigao Road, Hi-Tech Zone, Weihai, Shandong, China

TEL: +86 631 5660699 FAX: +86 631 5622563 P.C.: 264210

E-mail: info@weigaogroup.com exp@weigaogroup.com

Http: [//www.wegomedical.com/](http://www.wegomedical.com/)

This copyright belongs to Wego Medical International Co., Ltd. Our company reserves the right to investigate legal liability in case of infringement.

2024.2 EDITION

CE 0123

FDA





Sensory Control Protection



WG07-001 Surgical gown

- Detail information of exhibits
 - The clothes worn by surgeons and other medical staff during the operation in the operating room. Disposable surgical gowns are used to protect patients from infections by medical staff during the operation; they can also effectively protect medical staff.
 - Surgical gown is sewed by ultrasonic technology, without needle eyes, to ensure the antibacterial and anti-permeability of the surgical gown.
 - Reinforced surgical gown adds a chest sticker and two sleeve stickers on the basis of the standard surgical gown, which enhances the barrier performance of the surgical gown (high-risk parts) against bacteria and liquids.
 - Meet the requirements of EU standards (EN13795), standard (YY/T0506) standards and other performance.

WG07-002 Medical Cap

- Material: PP, SMS Non Woven
- Color: Blue, Green, Pink, White, Black, Red etc.
- Size: 18" 21" 24" / customization



WG07-003 Medical Face Mask

- Prevents direct transmission or spread of pathogenic microorganisms, body fluids or particulate matters. Provides infection prevention and outstanding protection

- YY 0469-2011 Surgical mask
- YY/T 0969-2013 Single-use medical face mask
- GB/T 32610-2016 Daily protective mask



WG07-004 KF94 Mask

- High protection level, good air tightness, with breathing valve
- Size: 17.5x9.5cm
- Safety standard: ISO 13485
- Filter Rating: Grade two
- Color: White/Blue





WG07-005 Medical Isolation Mask

- Outstanding protection with splash-resistance of body fluids or blood
- Consists of protective cover, foam strip and fixing device. The protective cover is made of polycarbonate, the foam strip is made of polyurethane, and the fixing device is composed of a headband made of polyester or nylon, latex wire, and a square buckle made of polypropylene. Non-sterile, one-time use.
- In the medical institution, it plays a protective role during inspection and treatment, blocking body fluids, blood splashes or splashes.
- PET protective sheet, high transparency and clear vision.
- The mask is made of anti-fog material, which can effectively prevent the formation of water mist.
- Using ergonomic design, lined with polyurethane foam strips, soft and comfortable to wear.
- It is convenient and firm to wear, has a large protection area, and the protection range can reach the entire face.



WG07-006 Medical Isolation Eye Mask

- Outstanding protection with splash-resistance of body fluids or blood
- It is composed of protective cover, lens and headband. The protective cover is made of polyurethane and polypropylene, the lens is made of polycarbonate, and the headband is made of polyester nylon, latex wire.
- Non-sterile, one-time use.
- In the medical institution, it plays a protective role in checking and treating, blocking body fluids, blood splashes or splashes.



WG07-007 Isolation Gown

- Infection prevention and outstanding protection
- Using non-woven fabric or film-coated non-woven fabric as the main raw material, made by cutting and sewing. Non-sterile, single use. It is used for general isolation in outpatient clinics, wards, inspection rooms, etc. of medical institutions.
- Multiple protection, light and comfortable, breathable and soft.
- Has good water resistance.
- Anti-fine particles, effective isolation.
- The head is designed with a cap and the cap is elastic to effectively block bacteria and particles.
- One-piece zipper design, easy to put on and take off, and better prevent bacteria and particles from entering.
- The cuffs are designed with elasticated mouth, which can be stretched freely.
- Elastic waist design to fit the body to meet the needs of different body shapes.



WG07-008 Medical Protective Gown

- Infection prevention and outstanding protection
- It is made of non-woven material with breathable film, the structure is one-piece, consisting of a hooded top, pants, and foot covers (optional), the cuffs, ankle openings, and foot cover openings are elastic closures, and the hat face closure and waist closure Use elastic closing. It is used to provide barrier and protection for the potentially infectious patients' blood, body fluids, secretions, etc. that clinical medical personnel come into contact with at work.
- In line with the national standard GB19082-2009 "Technical Requirements for Medical Disposable Protective Clothing", it is antibacterial, antistatic, and effective in blocking liquids and particles.
- Soft, light, breathable and comfortable to wear.
- The surface is smooth and not easy to get wet, which can effectively prevent the adhesion of bacteria.