

WEGO

LEADING BRAND OF MEDICAL
DEVICES IN CHINA



WEIGAO MEDICAL INTERNATIONAL CO., LTD.

ADD: No.1 Weigao Road, Hi-Tech Zone, Weihai, Shandong, China

TEL: +86 631 5660699 FAX: +86 631 5622563 P.C.: 264210

E-mail: info@weigaogroup.com exp@weigaogroup.com

Http: [//www.wegomedical.com/](http://www.wegomedical.com/)

This copyright belongs to Wego Medical International Co., Ltd. Our company reserves the right to investigate legal liability in case of infringement.

2024.2 EDITION

CE 0123

FDA





Medication Injection

Disposable Syringe (WG01-003-001)

- Material: Polypropylene
- Parts: Barrel+Piston+Plunger+Needle(or w/o needle) Luer Lock or Luer Slip
- Packing: Blister/PE
- Piston materials: Latex, Polyisoprene, Butyl
- Sterilization: EO

Code Number	Volume (ml/cc)	Needle		Luer Lock	Luer Slip	Packing
WGS1WL	1	W		L		Blister/ PE
WGS1WOS			W/O		S	Blister/ PE
WGS1WS		W			S	Blister/ PE
WGS1WOL			W/O	L		Blister/ PE
WGS2WL	2	W		L		Blister/ PE
WGS2WOS			W/O		S	Blister/ PE
WGS2WS		W			S	Blister/ PE
WGS2WOL			W/O	L		Blister/ PE
WGS3WL	3	W		L		Blister/ PE
WGS3WOS			W/O		S	Blister/ PE
WGS3WS		W			S	Blister/ PE
WGS3WOL			W/O	L		Blister/ PE
WGS5WL	5	W		L		Blister/ PE
WGS5WOS			W/O		S	Blister/ PE
WGS5WS		W			S	Blister/ PE
WGS5WOL			W/O	L		Blister/ PE
WGS10WL	10	W		L		Blister/ PE
WGS10WOS			W/O		S	Blister/ PE
WGS10WS		W			S	Blister/ PE
WGS10WOL			W/O	L		Blister/ PE
WGS20WL	20	W		L		Blister/ PE
WGS20WOS			W/O		S	Blister/ PE
WGS20WS		W			S	Blister/ PE
WGS20WOL			W/O	L		Blister/ PE
WGS30WL	30	W		L		Blister/ PE
WGS30WOS			W/O		S	Blister/ PE
WGS30WS		W			S	Blister/ PE
WGS30WOL			W/O	L		Blister/ PE
WGS50WL	50	W		L		Blister/ PE
WGS50WOS			W/O		S	Blister/ PE
WGS50WS		W			S	Blister/ PE
WGS50WOL			W/O	L		Blister/ PE
WGS60WL	60	W		L		Blister/ PE
WGS60WOS			W/O		S	Blister/ PE
WGS60WS		W			S	Blister/ PE
WGS60WOL			W/O	L		Blister/ PE
WGS100WL	100	W		L		Blister/ PE
WGS100WOS			W/O		S	Blister/ PE
WGS100WS		W			S	Blister/ PE
WGS100WOL			W/O	L		Blister/ PE



Disposable Insulin Syringe (WG01-003-002)

- Material: Polypropylene
- Parts: Barrel+Piston+Plunger+Fixed Needle+Cap
- Packing: Individual Blister/Self-Sealed Package
- Piston materials: Latex, Polyisoprene, Butyl
- Sterilization: EO

Code Number	Volume (ml/cc)	Needle Size		Packing	Colour Identification
		OD	L		
WGIN5J	1	29G	1/2	Blister/ PE	Orange
WGIN5S	0.5	30G	3/8		
WGIN5J	0.3	31G	5/16		
WGIN3J	0.3	32G	5/16	Blister/ PE	Red
WGIN1J	1	29G	1/2		
WGIN5J	0.5	30G	3/8		
		31G	5/16		
		32G	5/16		

* Customized needle size is available





Disposable Perfusor Syringe (WG01-003-003)

- Material: Polypropylene
- Parts: Barrel+Piston+Plunger
- Packing: Blister/PE
- Piston materials: Latex, Polyisoprene, Butyl
- Sterilization: EO
- Type: Conventional/Light-Resistant
- Pressure: 200pa
- Size: 50/60 ML(cc)

Disposable Pump Syringe (WG01-003-007)

- Material: Polypropylene
- Parts: Barrel+Piston+Plunger
- Packing: Blister/PE
- Piston materials: Latex, Polyisoprene, Butyl
- Sterilization: EO

Code Number	Volume (ml/cc)	Needle		Luer Lock	Packing
WGPS3	3	W	W/O	L	Blister/ PE
WGPS5	5	W	W/O	L	Blister/ PE
WGPS10	10	W	W/O	L	Blister/ PE
WGPS20	20	W	W/O	L	Blister/ PE
WGPS30	30	W	W/O	L	Blister/ PE
WGPS50	50	W	W/O	L	Blister/ PE
WGPS60	60	W	W/O	L	Blister/ PE



Pre-filled Flush Syringe (WG01-003-004)

- Sterilization method: high temperature moist heat sterilization
- Main materials: polypropylene, butyl rubber

Specification	Single PE bag	Box size	Capacity of a box	Carton size	G.W
3ml	30pcs/box	26.6*19*9.1cm	480pcs/carton	530*385*390mm	8kg
5ml	30pcs/box	22.5*16*9.5cm	480pcs/carton	475*345*415mm	10kg
10ml	30pcs/box	22*22*9.5cm	480pcs/carton	475*465*415mm	14kg



Disposable Sterile Rinser with Catheter Tip (WG01-003-005)

- Material: Polypropylene
- Parts: Barrel+Piston+Plunger+Cap
- Packing: Blister/PE
- Piston materials: Latex, Polyisoprene, Butyl
- Sterilization: EO
- Size: 50/60 ML(cc)



Disposable Light-resistant Syringe (WG01-003-008)

- Material: Polypropylene
- Parts: Barrel+Piston+Plunger
- Packing: Blister/PE
- Piston materials: Latex, Polyisoprene, Butyl
- Sterilization: EO

Code Number	Volume (ml/cc)	Needle		Luer Lock	Packing
WGLRS3	3	W	W/O	L	Blister/ PE
WGLRS5	5	W	W/O	L	Blister/ PE
WGLRS10	10	W	W/O	L	Blister/ PE
WGLRS20	20	W	W/O	L	Blister/ PE
WGLRS30	30	W	W/O	L	Blister/ PE
WGLRS50	50	W	W/O	L	Blister/ PE

Feeding Syringe (WG01-003-006)

Feeding syringe are designed specifically for enteral processes. These processes include initial tube placement, flushing, irrigation and more. The ENFit connector reduces the risk of misconnections to the tubing. Also, the body is clear for easy measuring against the clearly marked graduated length markings.

- Material: Polypropylene
- Parts: Barrel+Piston+Plunger+Cap
- Packing: Blister/PE
- Piston materials: Latex, Polyisoprene, Butyl
- Sterilization: EO
- Size: 3/5/10/20/30/50 ML(cc)



Disposable Fixed Dose Vaccine Syringe (WG01-003-009)

- Material: Polypropylene
- Parts: Barrel+Piston+Plunger+Fastener+Fixed Needle+Cap
- Packing: Blister/PE
- Piston materials: Latex, Polyisoprene, Butyl
- Sterilization: EO

Code Number	Volume (ml/cc)	Needle Size		Packing
		OD	L	
WGV501	1	26G	1/2"	Blister/ PE
WGV502	0.9			Blister/ PE
WGV503	0.8			Blister/ PE
WGV504	0.7			Blister/ PE
WGV505	0.6	25G	5/8"	Blister/ PE
WGV506	0.5			Blister/ PE
WGV507	0.4	24G	1"	Blister/ PE
WGV508	0.3			Blister/ PE
WGV509	0.2	23G	1"	Blister/ PE
WGV51	0.1			Blister/ PE

