

# WEGO

LEADING BRAND OF MEDICAL  
DEVICES IN CHINA



**WEIGAO MEDICAL INTERNATIONAL CO., LTD.**

ADD: No.1 Weigao Road, Hi-Tech Zone, Weihai, Shandong, China

TEL: +86 631 5660699 FAX: +86 631 5622563 P.C.: 264210

E-mail: [info@weigaogroup.com](mailto:info@weigaogroup.com) [exp@weigaogroup.com](mailto:exp@weigaogroup.com)

Http: [//www.wegomedical.com/](http://www.wegomedical.com/)

This copyright belongs to Wego Medical International Co., Ltd. Our company reserves the right to investigate legal liability in case of infringement.

2024.2 EDITION





Surgical Supply



#### WG05-001-001 Disposable Surgical Blade

Surgical Blades are critically made in accordance with ISO9001/ISO7740 standard.

They supply very high quality to high level doctors worldwide. Our Surgical Blades have the most popular sizes to satisfy different surgery demands

Material: carbon steel or stainless steel


Size: 10#-36#



#### WG05-001-002 Disposable Skin Graft Blade

Mdterial: carbon steel or stainless steel

Size: 160mm



#### WG05-001-003 Disposable Stitch Cutter Blade

Special Blades include the following sizes:stitch cutter blade 65# and 110#, These blades can meet your different surgical medical use requirements.

Raw material:carbon steel or stainless steel

Size: 65mm, 110mm



#### WG05-001-004 Disposable Medical Razor

Carbon steel or stainless steel,single or double blade.

All razors come with the highest quality Stainless Steel blades or carbon steel,Comfart-Coated for smooth shaves.



#### WG05-001-005 Disposable Surgical Scalpel

Material:Stainless Steel

Size:Customize

Color:Silver,Golden,Black,Blue Coated,Multi Color,Color Coated Or Customized



Finish: Mirror,Dull,Satin

Design: Custom Design



Sterile Surgical Suture

Wegosutures covers all the suture lines including Absorbable and Non-Absorbable sutures approved by ISO, CE and FDA, marketing globally under the umbrella of Product Liability Insurance since 2007. Strong and sharp needle combined well with soft and smooth thread with sophisticated manufacturing process, offering various easy-to-use packaging that makes Wegosutures become the reliable choice of surgery.





Surgical Needle

Medical grade alloy was made into atraumatic suture needle by the machines and technology originate from US and Europe. Precision drillings always centralize with wire brings lower spoil rate during the crimping procedure. Carefully designed production procedures well balanced the Ductility and Bend Moment with the excellent Penetration. Medical silicone coating ensure the smooth penetration of several pass.



Category	Alloy	Needle Tip	Curve circle	length	Bore Range	Hole Range
Drilled-End	420	TA TC BL RC CU P	1/4,1/2,3/8,5/8,Straight,Ski,J(Hook),Compound	10-300mm	0.33-1.6mm	0.1-1.2mm
	455	TA TC BL RC CU P	1/4,1/2,3/8,5/8,Straight,Ski,J(Hook),Compound	10-90mm	0.33-1.3mm	0.2-1.05mm
	470	TA TC BL RC CU P	1/4,1/2,3/8,5/8,Straight,Ski,J(Hook),Compound	10-90mm	0.33-1.3mm	0.15-0.85mm
	302/304	TA TC BL RC CU P	1/4,1/2,3/8,5/8,Straight,Ski,J(Hook),Compound	16-48mm	0.41-1.7mm	0.2-1.5mm
Regular Eyed	420	TA CU	1/2,3/8,Straight	10-90mm	0.2-1.3mm	
	304	TA CU	1/2,3/8,Straight	10-90mm	0.2-1.3mm	

TA- Taper Point / TC- Taper Cut / BL- Blunt / RC- Reverse Cutting / CU- Conventional Cutting /Cutting / P- Extra Sharp  
\* - contact us for detailed available specification

Cassette Suture

Medical grade compound becomes strong and clear cassette after precision mould injection.15-100 meters sutures in the cassette cut off the cost of bulk surgeries in field. Ultrasonic Sealing protects the cassette out of leakage. Well fit with the racks in the market easy to carry.



Product	absorption	Color	USP Range
WEGO PGA Cassette	Synthetic Absorbable	Violet	5#-3/0
WEGO PDO Cassette	Synthetic Absorbable	Violet	2#-3/0
WEGO Chromic Catgut Cassette	Absorbable	Brown	2#-3/0
WEGO Plain Catgut Cassette	Absorbable	Beige	2#-3/0
WEGO Nylon Cassette	Non-absorbable	Blue / Black	2#-3/0
WEGO SUPRAMID Cassette	Non-absorbable	Black/ White	7#-3/0
WEGO Silk Cassette	Non-absorbable	Balck	2#-3/0
WEGO Polyester Cassette	Non-absorbable	Green / White	5#-3/0

Non Sterile Surgical Suture Material

Available to supplying from Silk to PDO covering the most requirement of sutures manufacturing with USP 2# to 6/0. Unique Braided Structure on PGA thread offering comparable softness and strength performance with leading brands. We offer special customized service on both Braided and Monofilament sutures.



Structure	Product	absorption	Material	Coating	Color	USP Range
Braided	PGA	Synthetic Absorbable	Polyglycolic Acid	Ca-Sterate	Violet / Undyed	5#-6/0
Monofilament	PDO	Synthetic Absorbable	Polydioxanone	N/A	Violet	3#-6/0
Monofilament	PGCL	Synthetic Absorbable	Poly(glycolide-co-caprolactone)(75/25)	N/A	Undyed / Violet	2#-6/0
Monofilament	Polypropylene	Non-absorbable	Polypropylene	N/A	Blue / Clear	3#-6/0
Monofilament	Nylon	Non-absorbable	Polyamide 6 and 66	N/A	Blue / Black /Clear	7#-6/0
Braided	Silk	Non-absorbable	100% Protein fiber spun by silkworm Bombyx Mori	Silicone	Balck	2#-6/0
Braided	Nylon	Non-absorbable	Polyamide 66	Uncoated Silicone	White/Balck	2#-6/0
Braided	Polyester	Non-absorbable	Polyethylene terephthalate	Silicone	Green	2#-6/0

Medical Grade Chemical Compound

Wego PVC compound was made by top grade resin and materials in clean room. Strictly quality assurance system brings low rate impurities and high consistence on crystallization during injection and extrusion process, this gains the medical device better stability in shelf life and performance.



Material	Type	Code	Main Application
Medical Grade PVC Compound	MT	T0C IA	Extrusion Catheter/Tube
	MD	D0I IA	Injection Drip Chamber
	MT	P0I IA	Injection Syringe Gasket
	MT	T0LEIA	Extrusion Elastic Catheter/Tube
	MT	T0HEIA	Extrusion High Elastic Catheter/Tube/Capillary Catheter
	MD	DOB IA	Extrusion Blister Drip Chamber
	MD	F0I ID	Filters and Connectors by Injection
	MD	N0I ID	Accessories of medical device by injection
	MT	TIC IA	DINP Plasticizer, Extrusion Catheter/Tube
	MD	DII IA	DINP Plasticizer, Injection Drip Chamber
Medical Grade TPE Granula	MT	T9E IA-H	Extrude the TPE Capillary Catheter
	MT	T9E IA	Extrude the TPE Catheter/Tube
	MD	D9I ID	Injection TPE Drip Chamber